

L-2000™

Drinking Valve



Eliminate Manual Delivery of Drinking Water

Today's ventilated caging systems provide maximum housing density and improved environmental control to extend the length of time between cage changes. The labor saving advantages of these systems are drastically reduced when water bottles are used. Bottles require frequent changing, filling, washing, and material handling.

In order to take full advantage of these advanced housing systems, the reliable L-2000™ bedding resistant drinking valve, combined with our ReCIRC™ water distribution system streamlines workflow while providing a continuous supply of fresh contaminant free drinking water. The activation force required to operate this valve is designed for adult and weanling rodents. This valve has proven its reliability for over 14 years. We supply both cage mounted valves or rack mounted valves with easy disconnect for cleaning.

Precision manufactured of 316 stainless steel with an operating range of 0.5 to 5 psi. The L-2000™ valve and Slim Line™ stainless steel rack manifolds are available on all cage manufacturers systems. Be sure to specify our system when ordering your animal caging systems. Our animal watering solutions also include ReCIRC™ Room Distribution Systems, providing continuous circulation of contaminant free water, using UV sterilization and RO filtration.



SYSTEMS ENGINEERING
SE LAB GROUP INC.

932 Kaiser Road, Napa, CA 94558 • 877-253-8852 • www.selabgroup.com
© SE Lab Group Inc. 2003 ReCIRC, ReCIRC Rack are trademarks of SE Lab Group Inc.